

# **REQUISITION FORM** Telomere Length Measurements

Today's date:			Store patient sample at room temperature  Do not refrigerate						
PATIENT INFORMATION									
Patient's last name:	ast name: First: Middle:			Birth Date: mm / dd / yyyyy		Sex □ M □ F			
Patient ID#:				Sample Collection Date		Time hh / mm			
REASON FOR TESTING									
☐ Bone Marrow Failure	☐ Imr	nunodeficiency	☐ Lymphoid Ma	hoid Malignancy					
☐ Pulmonary Fibrosis	☐ Oth	ner Lung Disease	☐ Other, please	ease specify:					
ORDERING INFORMATION									
Physician:				Dept.:					
Hospital:									
Address:									
City:				Prov:	Pos	stal Code:			
The person listed as the Ordering Physician is authorized by law to order the test.				Results to be sent by:					
Authorized Signature (Req		☐ Fax:							
		TEC	T DEQUESTS			_			
TEST REQUESTED  Repeat Diagnostics uses the Flow FISH procedure.  Turnaround time is within 3 weeks. For expedite service, please contact us.									
_'		th measurements for total l <b>ym</b>			•	• •			
l <u> </u>	_				-	T-cells and NK cells. \$1050.00			
☐ Medical Consultation		0 for a written evaluation by a							
information, such as family history, clinical history, current working diagnosis, symptoms and lab investigations. If the space allocated is not enough, please provide additional information on a separate sheet:									
PATIENT MEDICAL INFORMATION									
BILLING OPTIONS									
	Institu	(We do not invoice utional Billing:	neamicare insura	· ,	ng Credit	card (VISA & MasterCard)			
Hospital:		<b>.</b>		Name on Credit Card:					
Department:				Address:					
Contact:				City:					
Address:				Prov:	Pos	stal Code:			
City:				Card number:					
Prov:		Postal Code:		Exp. Date (mmyy):		CVC:			
Tel:				Signature of Cardholder:					
Email:			Please charge the above credit card in the amount of \$						



# TELOMERE LENGTH MEASUREMENTS SPECIMEN COLLECTION AND SHIPPING PROCEDURE

#### **BEFORE COLLECTION OF BLOOD**

Sample should <b>only</b> be collected and shipped on Monday, Tuesday or Wednesday.	
Requisition Form check list	
☐ Patient name is filled in and matches blood tube ID (first identifier)	
☐ Second patient identifier (date of birth, unique ID number) is filled in and matches blood tu	be
☐ Ordering information is complete and signed by the requesting physician	
Result send out information is completed	
Assay type (2 or 6-panel) and optional consultation are selected accordingly	
☐ Payment information is completed	
SPECIMEN COLLECTION	
□ Label the specimen tube with: □ Patient Name and ID # □ Age □ Sex	15 to
□ Date and time of collection	N. 52 50,000 reserves

#### SPECIMEN PACKING AND SHIPPING

Collect blood in EDTA anti-coagulant tube.5-10ml of blood is required for successful testing.

## SHIPPING MATERIAL

SHIP	PING WATERIAL	
	UN3373 shipping box measuring approximately 9" X 4" X 4", labeled "Bi	ological Substance Category B)
	Specimen bag or sealable plastic bag.	
	Absorbent material such as paper towel.	-
	Packing tape.	OFFICIAL SUBSTANCE, CATEGORY &
	Address label.	TC-125-18
	FedEx Clinical Pak (provided free of charge from FedEx)	UN3973/
	Intra Canada air waybill.	100
	For more information on how to ship clinical samples visit FedEx at http://images.fedex.com/downloads/shared/packagingtips/pointers.pdf	

### **SHIPPING**

- 1. Place blood collection tube(s) in sealable plastic bag.
- 2. Place bag in shipping container. ICE PACKS ARE NOT REQUIRED

□ Store patient sample at **room temperature** until pick-up by courier.

□ All blood shipments to Repeat Diagnostics must arrive within 2 days and in good condition.

- 3. Place enough absorbent material in shipping container so that blood tubes do not roll around.
- 4. Seal shipping container with packing tape.
- 5. Attach address label to top of shipping container.
- 6. Place shipping container and requisition form inside FedEx Clinical Pak.
- 7. Fill out the Intra Canada air waybill form.
- 8. Ship on day of collection by FedEx Priority or FedEx First Overnight to:

Repeat Diagnostics Inc. Suite 309 - 267 West Esplanade North Vancouver, BC V7M 1A5 Canada

9. Inform Repeat Diagnostics by email at test@repeatdx.com of date shipped and tracking number.