

Today's date:		Store patient sample at room temperature Do not refrigerate	
PATIENT INFORMATION			
Patient's last name:		First:	Middle:
		Birth Date: mm / dd / yyyy	
		Sex: <input type="checkbox"/> M <input type="checkbox"/> F	
Patient ID#:		Sample Collection Date: mm / dd / yyyy	
		Time: hh / mm	
REASON FOR TESTING			
<input type="checkbox"/> Bone Marrow Failure	<input type="checkbox"/> Immunodeficiency	<input type="checkbox"/> Lymphoid Malignancy	<input type="checkbox"/> Myeloid Malignancy
<input type="checkbox"/> Pulmonary Fibrosis	<input type="checkbox"/> Other Lung Disease	<input type="checkbox"/> Other, please specify:	
ORDERING INFORMATION			
Physician:		Dept:	
Hospital:		Test PO Number:	
Address:			
City:		Country:	
The person listed as the Ordering Physician is authorized by law to order the test.		Results to be sent by:	
Authorized Signature (Required):		<input type="checkbox"/> Fax: _____	
		<input type="checkbox"/> Email: _____	
TEST REQUESTED			
Repeat Diagnostics uses the Flow FISH procedure. Turnaround time is within 3 weeks. For expedited service, please contact us.			
<input type="checkbox"/> 2-Panel Assay Telomere length measurements for total lymphocyte and granulocyte population only.			
<input type="checkbox"/> 6-Panel assay Telomere length measurements for total lymphocytes and granulocytes as well as B-cells , T-cells and NK cells .			
<input type="checkbox"/> Medical Consultation €250.00 for a written evaluation by a hematopathologist to accompany the test results. Provide pertinent patient information, such as family history, clinical history, current working diagnosis, symptoms and lab investigations. If the space allocated is not enough, please provide additional information on a separate sheet.			
PATIENT MEDICAL INFORMATION			
BILLING OPTIONS			
(We do not invoice healthcare insurance companies)			
Institutional Billing:		International Bank Transfer	
Hospital:		IBAN: NL77 INGB 0006 6060 80 / BIC: INGBNL2A	
Department:		Repeat Diagnostics Europe BV	
Contact:		Heresingel 4	
Address:		9711 ES Groningen	
City:		The Netherlands	
Country:	Postal Code:		
Tel:			
Email:			

BEFORE COLLECTION OF BLOOD

Sample should **only** be collected and shipped on Monday, Tuesday or Wednesday.

Requisition Form check list

- Patient name is filled in and matches blood tube ID (first identifier).
- Second patient identifier (date of birth, unique ID number) is filled in and matches blood tube.
- Ordering information is complete and signed by the requesting physician.
- Result send out information is completed.
- Assay type (2 or 6-panel) and optional consultation are selected accordingly.
- Payment information is completed.

SPECIMEN COLLECTION

- Label the specimen tube with:
 - Patient Name and ID #
 - Age
 - Sex
 - Date and time of collection
- Collect blood in EDTA anti-coagulant tube.
- 5-10ml of blood is required for successful testing.
- Store patient sample at room temperature until pick-up by courier.
- All blood shipments to Repeat Diagnostics must arrive within 2 days and in good condition.

**SPECIMEN PACKING AND SHIPPING**SHIPPING MATERIAL

- UN3373 shipping box measuring approximately 9" X 4" X 4", labeled "Biological Substance Category B".
- Specimen bag or sealable plastic bag.
- Absorbent material such as paper towel.
- Packing tape.
- Address label.

SHIPPING

1. Place blood collection tube(s) in sealable plastic bag.
2. Place bag in shipping container. ICE PACKS ARE NOT REQUIRED.
3. Place enough absorbent material in shipping container so that blood tubes do not roll around.
4. Seal shipping container with packing tape.
5. Attach address label to top of shipping container.
 - If necessary, add documents for customs. The required EORI number is 5379687.
 - Please write the following sentence clearly on the package: "Eilige Patientenprobe, bitte sofort ausliefern"
6. Ship samples by express delivery on day of collection to:

Telomerdiagnostik
Frau Lucia Vankann
Hämatologisches Labor
Klinik für Hämatologie, Onkologie, Hämostaseologie und Stammzelltransplantation
Pauwelsstrasse 30
52074 Aachen

7. Inform Repeat Diagnostics by email at request@repeatdx.com of date shipped and tracking number.